
CLEARVIEW GLASS CALIBRATION CYLINDER CAD FILE CODING

FILE CODES Example: CLEARVIEW GLASS CYLINDER, 2000 mL, PVC INSERT WITH 1" THREADED CONNECTION



CYLINDER SIZE

100 = 100 mL CYLINDER

250 = 250 mL CYLINDER

500 = 500 mL CYLINDER

1000 = 1000 mL CYLINDER

2000 = 2000 mL CYLINDER

4000 = 4000 mL CYLINDER

INSERT MATERIAL

PV = PVC / CPVC (CORZAN)

POLYVINYLCHLORIDE /

CHLORINATED

POLYVINYLCHLORIDE

= PP / PVDF

POLYPROPYLENE /

POLYVINYLIDENE FLUORIDE

SS = 316 STAINLESS STEEL

CONNECTION SIZE

50 = 1/2" (DN 15)

STANDARD ON 100 mL TO 1000 mL SIZES

75 = 3/4" (DN 20)

OPTION FOR 500 & 1000 mL SIZES

10 = 1" (DN 25)

STANDARD ON 2000 mL & 4000 mL SIZES

OPTION FOR 500 mL & 1000 mL SIZES

15 = 1-1/2" (DN 40)

OPTION FOR 2000 mL & 4000 mL SIZES

20 = 2" (DN 50)

OPTION FOR 4000 mL SIZE

CONNECTION TYPE

N = THREADED /

SOCKET

F = FLANGED

DO NOT COPY WITHOUT PRIOR WRITTEN APPROVAL FROM PRIMARY FLUID SYSTEMS INC.

PRIMARY FLUID SYSTEMS INC.

1050 Cooke Blvd., Burlington, Ontario L7T 4A8

CLEARVIEW GLASS CYLINDERS

DATE DRAWN: 1/23/2018 SCALE: N/A SHEET: 1 DF 1

PROJECT: CVG DRAWN BY: F:

0 CODING FOR WEBSITE -- -REV DESCRIPTION DATE NAME

NOTE: OVERALL WIDTH, HEIGHT, AND DEPTH OF DRAWING TOLERANCE: ± 1/4"

APPROVED CHECKED DWG NO.:
BY: CVG-CAD-FILE-CODING 0